

Bioinformatics Lessons Schedule

- RNA-seq
- single cell RNA-seq
- RRBS

| Date | Subject |
|-------|----------------------------|
| 12-24 | Christmas break |
| 12-31 | Christmas break |
| 01-07 | Process RNA-seq |
| 01-14 | Process RNA-seq, continued |
| 01-21 | Process RNA-seq, continued |
| 01-28 | Analyze RNA-seq |
| 01-28 | Analyze RNA-seq, continued |

RNA-seq

Quality Check

FastQC

- Before going forward, we want to check the quality of the data
 - How much did the sequencer fail?
 - Did we sequence mostly our sample DNA?
- FastQC is a program from the Babraham Institute in the UK that creates an html report on the quality of the sequencing data
 - Has 11 quality control checks that it does

Basic Statistics

Good Quality



Basic Statistics

| Measure | Value |
|-----------------------------------|-------------------------|
| Filename | good_sequence_short.txt |
| File type | Conventional base calls |
| Encoding | Illumina 1.5 |
| Total Sequences | 250000 |
| Sequences flagged as poor quality | 0 |
| Sequence length | 40 |
| %GC | 45 |

Bad Quality




Basic Statistics

| Measure | Value |
|-----------------------------------|-------------------------|
| Filename | bad_sequence.txt |
| File type | Conventional base calls |
| Encoding | Illumina 1.5 |
| Total Sequences | 395288 |
| Sequences flagged as poor quality | 0 |
| Sequence length | 40 |
| %GC | 47 |

Per base sequence quality


Good Quality

✔ Per base sequence quality



Bad Quality


✖ Per base sequence quality



Per tile sequence quality


Good Quality

✔ Per tile sequence quality



Bad Quality


✘ Per tile sequence quality



Per sequence quality scores


Good Quality

✓ Per sequence quality scores



Bad Quality


✓ Per sequence quality scores



Per base sequence content


Good Quality

✔ Per base sequence content



Bad Quality


⚠ Per base sequence content



Per sequence GC content


Good Quality

✓ Per sequence GC content



Bad Quality


! Per sequence GC content



Per base sequence quality


Good Quality

✓ Per base N content



Bad Quality


✓ Per base N content



Per base sequence quality


Good Quality

✔ Sequence Length Distribution



Bad Quality


✔ Sequence Length Distribution



Sequence Duplication Levels


Good Quality

✔ Sequence Duplication Levels



Bad Quality

! Sequence Duplication Levels



Overrepresented sequences

Good Quality



Overrepresented sequences

No overrepresented sequences


Bad Quality

| Overrepresented sequences | | | | |
|---|-------|----------------------|--|-----|
| Sequence | Count | Percentage | Possible Source | |
| AGAGTTTATCGCTTCACGAGCGCAGAGTTTACACCTTC | 2065 | 0.5224039181558763 | No | Hit |
| GATGCGGCTATCAACTCGAGGTTTATCGCTTCATG | 2047 | 0.5178592762542754 | No | Hit |
| ATTGGGCTATCAACTCGAGGTTTATCGCTTCATG | 2014 | 0.5095019327680071 | No | Hit |
| CGATAAAATGATGGGCTATCAACTCGAGGTTTAT | 1913 | 0.4839509420979134 | No | Hit |
| GTATCAACTCGAGGTTTATCGCTTCATGAGCGAGA | 1879 | 0.4753496185060066 | No | Hit |
| AAAAATGATGGGCTATCAACTCGAGGTTTATCGCT | 1846 | 0.4670012750197325 | No | Hit |
| TGATGGGCTATCAACTCGAGGTTTATCGCTTCAT | 1841 | 0.46573637449150995 | No | Hit |
| AACCTCGAGGTTTATCGCTTCATGAGCGAGTTAA | 1836 | 0.46447147396328753 | No | Hit |
| GATAAAATGATGGGCTATCAACTCGAGGTTTATC | 1831 | 0.4632065734350651 | No | Hit |
| AAATGATGGGCTATCAACTCGAGGTTTATCGCTTC | 1779 | 0.45005160794155147 | No | Hit |
| ATGATGGGCTATCAACTCGAGGTTTATCGCTTCA | 1779 | 0.45005160794155147 | No | Hit |
| AAATGATGGGCTATCAACTCGAGGTTTATCGCTTC | 1760 | 0.4452449859343061 | No | Hit |
| AAAAATGATGGGCTATCAACTCGAGGTTTATCGCT | 1729 | 0.4374026026593269 | No | Hit |
| CGTATCAACTCGAGGTTTATCGCTTCATGAGCGAGA | 1713 | 0.4333549299691494 | No | Hit |
| ATCCAACTCGAGGTTTATCGCTTCATGAGCGAGAG | 1708 | 0.4320902044079253 | No | Hit |
| CAGAGTTTATCGCTTCATGAGCGAGGTTTACACTT | 1684 | 0.42601849790532476 | No | Hit |
| TGACAGGTTTATCGCTTCATGAGCGAGGTTTAACT | 1668 | 0.4219708162150128 | No | Hit |
| CAACTCGAGGTTTATCGCTTCATGAGCGAGGTTA | 1668 | 0.4219708162150128 | No | Hit |
| TATCCAACTCGAGGTTTATCGCTTCATGAGCGAGA | 1630 | 0.4123575722005221 | No | Hit |
| GTATGGAAGCGATAAACTCTGAGGTTGATACGCCA | 1620 | 0.40982777114407726 | No | Hit |
| AACCTCTGAGGTTTATCGCTTCATGAGCGAGTTGG | 1616 | 0.4088158507214993 | No | Hit |
| CGAGGTTTATCGCTTCATGAGCGAGGTTTAACTT | 1580 | 0.39970856691829754 | No | Hit |
| TGGGCTATCAACTCGAGGTTTATCGCTTCATGAGC | 1569 | 0.3969257857562082 | No | Hit |
| GGGCTATCAACTCGAGGTTTATCGCTTCATGAGC | 1542 | 0.39009532290380683 | No | Hit |
| ATAAAATGATGGGCTATCAACTCGAGGTTTATCG | 1481 | 0.37466353645949285 | No | Hit |
| ACCTCGAGGTTTATCGCTTCATGAGCGAGGTTAAC | 1479 | 0.37415757624820384 | No | Hit |
| ATGGAAGCGATAAACTCTGAGGTTGATACGCCA | 1452 | 0.3673271133958026 | No | Hit |
| GATAAACTCTGAGGTTGATACGCCAATCATTTTATC | 1420 | 0.35923175001517876 | No | Hit |
| CGTATGGAAGCGATAAACTCTGAGGTTGATACGCCA | 1412 | 0.3572079091700229 | No | Hit |
| ACTCTCGCTCATGGAAGCGATAAACTCTGAGGTTGA | 1368 | 0.344607678452166524 | No | Hit |
| TAACTCTCGCTCATGGAAGCGATAAACTCTGAGGTTG | 1363 | 0.34481188399344276 | No | Hit |
| CATGGAAGCGATAAACTCTGAGGTTGATACGCCA | 1333 | 0.337222480824108 | No | Hit |
| CGATAAACTCTGAGGTTGATACGCCAATCATTTTAT | 1304 | 0.32988605776041774 | No | Hit |
| TAAAAATGATGGGCTATCAACTCGAGGTTTATCGC | 1277 | 0.32305559490801644 | No | Hit |
| GGCTATCAACTCGAGGTTTATCGCTTCATGAGCGA | 1262 | 0.31926089332334906 | No | Hit |
| TGGCTATGGAAGCGATAAACTCTGAGGTTGATACGC | 1233 | 0.3119244702595688 | No | Hit |
| GGAGCGATAAACTCTGAGGTTGATACGCCAATCATT | 1182 | 0.2990224848717897 | No | Hit |
| AAGCGATAAACTCTGAGGTTGATACGCCAATCATTT | 1136 | 0.2873854000121431 | No | Hit |
| ACTCTGAGGTTGATACGCCAATCATTTTATGAGGCG | 1133 | 0.28646245496920956 | No | Hit |
| AAACTCTGAGGTTGATACGCCAATCATTTTATGAGG | 1131 | 0.28612049494939206 | No | Hit |
| AAACTCTGAGGTTGATACGCCAATCATTTTATGAGG | 1129 | 0.2856145392726316 | No | Hit |
| AAGCGATAAACTCTGAGGTTGATACGCCAATCATTT | 1113 | 0.2815668575823197 | No | Hit |
| ATAAACTCTGAGGTTGATACGCCAATCATTTTATCG | 1111 | 0.28106089737103074 | No | Hit |
| AACCTCGAGGTTGATACGCCAATCATTTTATGAGGCG | 1083 | 0.273977454612985 | No | Hit |
| CTCGAGGTTTATCGCTTCATGAGCGAGGTTTAACT | 1055 | 0.2668940114549392 | No | Hit |
| TCTCGGCTCATGGAAGCGATAAACTCTGAGGTTGATA | 947 | 0.23957216004533402 | No | Hit |
| TGGAAGCGATAAACTCTGAGGTTGATACGCCAATCAT | 946 | 0.23931917993968954 | No | Hit |
| TAAAACTCTGAGGTTGATACGCCAATCATTTTATGGA | 912 | 0.2307178563477768 | No | Hit |
| GAGCGATAAACTCTGAGGTTGATACGCCAATCATT | 888 | 0.224646333812309 | No | Hit |
| GGCTATGGAAGCGATAAACTCTGAGGTTGATACGCC | 805 | 0.20364898504381615 | No | Hit |
| GGCGATAAACTCTGAGGTTGATACGCCAATCATTTTAA | 785 | 0.19898938293992632 | No | Hit |
| TGCGGCTATCAACTCGAGGTTTATCGCTTCATGAGG | 784 | 0.198336409825818 | No | Hit |
| CTCTCGGCTATGGAAGCGATAAACTCTGAGGTTGAG | 762 | 0.192770840501103 | No | Hit |
| TGCAACTCTGAGGTTTATCGCTTCATGAGCGAGGTT | 752 | 0.1902410394445806 | No | Hit |
| CCAACTCTGAGGTTTATCGCTTCATGAGCGAGGTT | 744 | 0.18821719859950212 | No | Hit |
| TGATGGAAGCGATAAACTCTGAGGTTGATACGCCA | 665 | 0.16823177025358726 | No | Hit |
| TCTCGGCTCATGGAAGCGATAAACTCTGAGGTTGATAC | 627 | 0.15861852623909656 | No | Hit |
| CGTTCAGAGGTTTATCGCTTCATGAGCGAGGTTTAA | 613 | 0.15507680476007366 | No | Hit |
| CGGCTCAGAGGATATCGGAGGATCGGAAGGCGGCTCAGC | 599 | 0.15153508328105078 | Illumina Paired End PCR Primer 2 (96% over 25bp) | |
| TCTCGAGGTTGATACGCCAATCATTTTATGGAAGCGCG | 585 | 0.1479933618020279 | No | Hit |
| CGCTTAAAGCTACAGTTATATGCTCGGGGGGTTTTTTT | 552 | 0.13964501831575965 | No | Hit |
| CTCTCGAGGTTGATACGCCAATCATTTTATGGAAGCGCG | 532 | 0.1345854162028698 | No | Hit |
| TGGGCTATGGAAGCGATAAACTCTGAGGTTGATACGC | 515 | 0.13028475440691342 | No | Hit |
| TGCGAGGTTGATACGCCAATCATTTTATGGAAGCGCGC | 505 | 0.12775495335044852 | No | Hit |
| GCTTAAAGCTACAGTTATATGCTCGGGGGGTTTTTTT | 411 | 0.10397482341988626 | No | Hit |

Adapter Content


Good Quality

✔ Adapter Content



Bad Quality

✔ Adapter Content



Sidenote - Loops in Bash

What's a loop?

What's a loop?

A loop is a piece of code that repeatedly executes a command for a certain condition until that condition is no longer true.


What's a loop?

A loop is a piece of code that repeatedly executes a command for a certain condition until that condition is no longer true.

```
for i in *.fastq.gz;  
    do echo $i;  
done
```



What's a loop?

A loop is a piece of code that repeatedly executes a command for a certain condition until that condition is no longer true.

```
for i in *.fastq.gz; 
do echo $i;
done
```




What's a loop?


A loop is a piece of code that repeatedly executes a command for a certain condition until that condition is no longer true.

```
for i in *.fastq.gz; 
do echo $i; 
done
```

What's a loop?

A loop is a piece of code that repeatedly executes a command for a certain condition until that condition is no longer true.

```
for i in *.fastq.gz; 
do echo $i; 
done 
```



How do you write a loop in bash?

How do you write a loop in bash?

```
for i in *.fastq.gz;  
do echo $i;  
done
```


How do you write a loop in bash?

starting
a loop




```
for i in *.fastq.gz;  
do echo $i;  
done
```

How do you write a loop in bash?


starting
a loop

variable
name



```
for i in *.fastq.gz;  
do echo $i;  
done
```


How do you write a loop in bash?




The diagram illustrates the components of a bash for loop. Three red boxes are positioned above the code: 'starting a loop' points to 'for', 'variable name' points to 'i', and 'meeting this condition' points to the entire 'in *.fastq.gz;' line. A red bracket is placed under the 'in *.fastq.gz;' line.

```
for i in *.fastq.gz;  
do echo $i;  
done
```


How do you write a loop in bash?




How do you write a loop in bash?




How do you write a loop in bash?



How do you write a loop in bash?



How do you write a loop in bash?



Now go try
this loop on
the server

Run FastQC

1. Go to the RNA-seq data directory
2. Make a directory to put the FastQC reports into, `fastqc`
3. Run `fastqc` on the samples

```
for i in *.fastq.gz; do fastqc $i -o fastqc/;  
done
```

Trim Bad Quality Sequences

What is trimming and why do it?

What is trimming and why do it?

- Trimming removes sequencing adapters, bad quality sequences, and/or other biased sequence information

What is trimming and why do it?

- Trimming removes sequencing adapters, bad quality sequences, and/or other biased sequence information
- Why is that important?
 - Helps prevent incorrect base calls by removing poor quality information
 - Increases speed and accuracy of alignment by removing artificial sequences and low quality sequences

What is trimming and why do it?

- Trimming removes sequencing adapters, bad quality sequences, and/or other biased sequence information
- Why is that important?
 - Helps prevent incorrect base calls by removing poor quality information
 - Increases speed and accuracy of alignment by removing artificial sequences and low quality sequences
- Trimming does two complementary things:
 1. Removes any sequence information that comes from library preparation or sequencing
 2. Removes low quality bases / low quality reads

Trim Sequences

1. Go back up to the rnaseq directory
2. Make a folder to put the analysis results in, `analysis`
3. Make a folder inside the analysis folder to put the trimmed reads in, `analysis/01_trim`

Trim Sequences

```
for i in *R1.fastq.gz;
do trim_galore
    --paired
    --fastqc
    --illumina
    --output analysis/01_trim/
    --retain_unpaired
    -q 30
    $i
    ${i/R1/R2};
done
```


Trim Sequences

```
for i in *R1.fastq.gz;  
do trim_galore  
    --paired  
    --fastqc  
    --illumina  
    --output analysis/01_trim/  
    --retain_unpaired  
    -q 30  
    $i  
    ${i/R1/R2};  
  
done
```



loop condition

Trim Sequences

```
for i in *R1.fastq.gz;  
do trim_galore  
    --paired  
    --fastqc  
    --illumina  
    --output analysis/01_trim/  
    --retain_unpaired  
    -q 30  
    $i  
    ${i/R1/R2};  
  
done
```

loop condition

call the program

Trim Sequences

```
for i in *R1.fastq.gz;
```

loop condition

```
do trim_galore
```

call the program

```
--paired
```

reads are paired-end

```
--fastqc
```

```
--illumina
```

```
--output analysis/01_trim/
```

```
--retain_unpaired
```

```
-q 30
```

```
$i
```

```
${i/R1/R2};
```

```
done
```

Trim Sequences

```
for i in *R1.fastq.gz;
```

loop condition

```
do trim_galore
```

call the program

```
--paired
```

reads are paired-end

```
--fastqc
```

run FastQC again after trimming

```
--illumina
```

```
--output analysis/01_trim/
```

```
--retain_unpaired
```

```
-q 30
```

```
$i
```

```
${i/R1/R2};
```

```
done
```

Trim Sequences

```
for i in *R1.fastq.gz;
```

loop condition

```
do trim_galore
```

call the program

```
--paired
```

reads are paired-end

```
--fastqc
```

run FastQC again after trimming

```
--illumina
```

trim Illumina adapters

```
--output analysis/01_trim/
```

```
--retain_unpaired
```

```
-q 30
```

```
$i
```

```
${i/R1/R2};
```

```
done
```

Trim Sequences

```
for i in *R1.fastq.gz;
do trim_galore
  --paired
  --fastqc
  --illumina
  --output analysis/01_trim/
  --retain_unpaired
  -q 30
  $i
  ${i/R1/R2};
done
```

Annotations for the code:

- loop condition
- call the program
- reads are paired-end
- run FastQC again after trimming
- trim Illumina adapters
- output goes here

Trim Sequences

```
for i in $(ls *.fastq.gz);
```

loop condition

```
do trim_galore
```

call the program

```
--paired
```

reads are paired-end

```
--fastqc
```

run FastQC again after trimming

```
--illumina
```

trim Illumina adapters

```
--output analysis/01_trim/
```

output goes here

```
--retain_unpaired
```

keep reads where one mate fails
trimming but the other doesn't

```
-q 30
```

```
$i
```

```
${i/R1/R2};
```

```
done
```

Trim Sequences

```
for i in *R1.fastq.gz;
do trim_galore
  --paired
  --fastqc
  --illumina
  --output analysis/01_trim/
  --retain_unpaired
  -q 30
  $i
  ${i/R1/R2};
```

Annotations:

- loop condition
- call the program
- reads are paired-end
- run FastQC again after trimming
- trim Illumina adapters
- output goes here
- keep reads where one mate fails trimming but the other doesn't
- read files

done

Trim Sequences

```
for i in $(ls *.fastq.gz);
```

loop condition

```
do trim_galore
```

call the program

```
--paired
```

reads are paired-end

```
--fastqc
```

run FastQC again after trimming

```
--illumina
```

trim Illumina adapters

```
--output analysis/01_trim/
```

output goes here

```
--retain_unpaired
```

keep reads where one mate fails
trimming but the other doesn't

```
-q 30
```

Keep bases at this quality or above

```
$i
```

```
${i/R1/R2};
```

```
done
```

Trim Sequences

```
for i in $(ls *.fastq.gz);
```

loop condition

```
do trim_galore
```

call the program

```
--paired
```

reads are paired-end

```
--fastqc
```

run FastQC again after trimming

```
--illumina
```

trim Illumina adapters

```
--output analysis/01_trim/
```

output goes here

```
--retain_unpaired
```

keep reads where one mate fails
trimming but the other doesn't

```
-q 30
```

Keep bases at this quality or above

```
$i
```

```
${i/R1/R2};
```

By default bases
quality less than
20 will be
trimmed and if
the read falls
below 20 bp, it
will be discarded;
we set the
minimum quality
to be 30

```
done
```

Trim Sequences

```
for i in $(ls *.fastq.gz);
```

loop condition

```
do trim_galore
```

call the program

```
--paired
```

reads are paired-end

```
--fastqc
```

run FastQC again after trimming

```
--illumina
```

trim Illumina adapters

```
--output analysis/01_trim/
```

output goes here

```
--retain_unpaired
```

keep reads where one mate fails
trimming but the other doesn't

```
-q 30
```

Keep bases at this quality or above

```
$i
```

read files

```
${i/R1/R2};
```

By default bases
quality less than
20 will be
trimmed and if
the read falls
below 20 bp, it
will be discarded;
we set the
minimum quality
to be 30

```
done
```

Trim Command

```
for i in *R1.fastq.gz; do trim_galore  
--paired --fastqc --illumina --output analysis/01_trim/  
--retain_unpaired -q 30 $i ${i/R1/R2}; done
```